

Packaging



Polyethylene film with ethylene oxide indicator

BIOLENE® polyethylene packaging system has imprinted, all along the roll, a chemical indicator for ethylene oxide sterilization.

- It can be sealed according to the material's size.
- The sterile material can be stored for a long time.
- Cost-effective relation compared with others packaging systems.
- Rolls in a range from 4 cm to 50 cm width and 70 microns thick.

Pouch

- One side medical paper / one side transparent film (for ease of content and seal integrity identifying).
- Indicator imprints for steam and ethylene oxide sterilization.
- Rolls in a range from 5 cm to 30 cm width and 200 m large.

Nitrile Gloves for EtO Sterilization

Protect against chemical, acids, oils and detergents.

Telephone/fax:
(+54 11) 4308-4963
(+54 220) 4041-502

Address:
Saavedra y Ruta Prov. 21
(Alt. Km. 20,7 de ruta 21)
(B1761BWN) Pontevedra,
Buenos Aires, Argentina.

E-mail:
info@biolene.com
ventas@biolene.com

www.biolene.com

Quality Management System



CERTIFIED
according to
ISO 9001:2008

BIOLENE®
ETHYLENE OXIDE STERILIZATION